

## Neutral Atoms Image Generator

Generated by Doxygen 1.9.5



## Chapter 1

# README



## Chapter 2

t

<col style="width: 100%" > </colgroup> <tbody>

Neutral Atoms Image Generator Model for generating simulated images for neutral atoms in a grid, imaged by different cameras </tbody>

Loading...

Searching...

No Matches

Class Index

[C](#) | [E](#) | [I](#) | [T](#)

C

[Camera](#) (neutral\_atom\_imaging\_simulation.Camera)

[CMOSCamera](#) (neutral\_atom\_imaging\_simulation.Camera)

E

[EMCCDCamera](#) (neutral\_atom\_imaging\_simulation.Camera)


[Experiment](#) (neutral\_atom\_imaging\_simulation.Experiment)

I

[ImageGenerator](#) (neutral\_atom\_imaging\_simulation.ImageGenerator)

T

[TweezerArray](#) (neutral\_atom\_imaging\_simulation.Experiment)

Generated by  class="footer" width="104" height="31" alt="doxygen" /> 1.9.5



## Chapter 3

# Namespace Index

### 3.1 Namespace List

Here is a list of all documented namespaces with brief descriptions:

<a href="#">Camera</a>	Module that implements different cameras . . . . .	??
<a href="#">Experiment</a>	Module that implements different experiment setups . . . . .	??
<a href="#">ImageGenerator</a>	Module for generating images of neutral atoms in a grid . . . . .	??





## Chapter 4

# Hierarchical Index

### 4.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator . . . . .	??
ABC	
neutral_atom_imaging_simulation.Camera.Camera . . . . .	??
neutral_atom_imaging_simulation.Camera.CMOSCamera . . . . .	??
neutral_atom_imaging_simulation.Camera.EMCCDCamera . . . . .	??
neutral_atom_imaging_simulation.Experiment.Experiment . . . . .	??
neutral_atom_imaging_simulation.Experiment.TweezerArray . . . . .	??



## Chapter 5

# Class Index

### 5.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

<a href="#">neutral_atom_imaging_simulation.Camera.Camera</a>	Abstract camera class . . . . .	??
<a href="#">neutral_atom_imaging_simulation.Camera.CMOSCamera</a>	Use this camera if the generated images should look like they are taken by a CMOS camera .	??
<a href="#">neutral_atom_imaging_simulation.Camera.EMCCDCamera</a>	Use this camera if the generated images should look like they are taken by an electron multiplying charge-coupled device (EMCCD) camera . . . . .	??
<a href="#">neutral_atom_imaging_simulation.Experiment.Experiment</a>	Abstract experiment class . . . . .	??
<a href="#">neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator</a>	Main class for generating images . . . . .	??
<a href="#">neutral_atom_imaging_simulation.Experiment.TweezerArray</a>	Use this if a tweezer array is to be simulated . . . . .	??



## Chapter 6

# Namespace Documentation

### 6.1 Camera Namespace Reference

Module that implements different cameras.

#### 6.1.1 Detailed Description

Module that implements different cameras.

### 6.2 Experiment Namespace Reference

Module that implements different experiment setups.

#### 6.2.1 Detailed Description

Module that implements different experiment setups.

### 6.3 ImageGenerator Namespace Reference

Module for generating images of neutral atoms in a grid.

#### 6.3.1 Detailed Description

Module for generating images of neutral atoms in a grid.




## Chapter 7

# Class Documentation

### 7.1 neutral\_atom\_imaging\_simulation.Camera.Camera Class Reference

Abstract camera class.

Inheritance diagram for neutral\_atom\_imaging\_simulation.Camera.Camera:



#### Public Member Functions

- def [get\\_image\\_creation\\_method](#) (self)
- def [set\\_zernike\\_coefficients](#) (self, typing.Union[np.ndarray, typing.Tuple[int, int, int, int, int, int, int, int, int, int, int, int, int, int, int, int]] zernike\_coefficients)  
*Function for setting the zernike coefficients.*
- def [set\\_library](#) (self, ctypes.CDLL [library](#))  
*Function for setting the image generation library.*
- def [library](#) (self)  
*Function for getting the image generation library.*

#### Public Attributes

- **zernike\_coefficients**

#### 7.1.1 Detailed Description

Abstract camera class.

#### 7.1.2 Member Function Documentation

##### 7.1.2.1 get\_image\_creation\_method()

```
def neutral_atom_imaging_simulation.Camera.Camera.get_image_creation_method (
    self )
```

Reimplemented in [neutral\\_atom\\_imaging\\_simulation.Camera.EMCCDCamera](#), and [neutral\\_atom\\_imaging\\_simulation.Camera.CMOSCamera](#)

```
def neutral_atom_imaging_simulation.Camera.Camera.library (
    self )
```

## Returns

```
def neutral_atom_imaging_simulation.Camera.Camera.set_library (
    self,
    ctypes.CDLL library )
```

<i>library</i>	The image generation C library
----------------	--------------------------------

[illegible]



## Parameters

<i>zernike_coefficients</i>	<p>An array of exactly 15 values that represent the coefficients for</p> <ol style="list-style-type: none"> <li>1. piston</li> <li>2. y-tilt</li> <li>3. x-tilt</li> <li>4. defocus</li> <li>5. oblique astigmatism</li> <li>6. vertical astigmatism</li> <li>7. vertical coma</li> <li>8. horizontal coma</li> <li>9. vertical trefoil</li> <li>10. oblique trefoil</li> <li>11. primary spherical</li> <li>12. vertical secondary astigmatism</li> <li>13. oblique secondary astigmatism</li> <li>14. vertical quadrafoil</li> <li>15. oblique quadrafoil</li> </ol>
-----------------------------	--


The documentation for this class was generated from the following file:

- neutral\_atom\_imaging\_simulation/Camera.py

## 7.2 neutral\_atom\_imaging\_simulation.Camera.CMOSCamera Class Reference

Use this camera if the generated images should look like they are taken by a CMOS camera.

Inheritance diagram for neutral\_atom\_imaging\_simulation.Camera.CMOSCamera:



### Public Member Functions

- def `__init__` (self, typing.Tuple[int, int] resolution, float dark\_current\_sampling\_alpha=None, float dark\_current\_sampling\_beta=None, float quantum\_efficiency=None, float numerical\_aperture=None, float physical\_pixel\_size=None, float magnification=None, float bias\_clamp=None, float bias\_stdev=None, float row\_noise\_stdev=None, float column\_noise\_scale=None, float flicker\_noise\_scale=None, float preamp\_gain=None, float readout\_stdev=None, float exposure\_time=None, int binning=1)
- def `get_image_creation_method` (self)

*Function for acquiring the function handle of the library that is used to generate images using this camera.*

- `def apply_settings (self)`

*Function for relaying any settings changes to the library.*

## Public Attributes

- `dark_current_sampling_alpha`
- `dark_current_sampling_beta`
- `quantum_efficiency`
- `numerical_aperture`
- `physical_pixel_size`
- `magnification`
- `bias_clamp`
- `bias_stdev`
- `row_noise_stdev`
- `column_noise_scale`
- `flicker_noise_scale`
- `preampgain`
- `readout_stdev`
- `exposure_time`
- `binning`
- `resolution`
- `zernike_coefficients`

### 7.2.1 Detailed Description

Use this camera if the generated images should look like they are taken by a CMOS camera.

### 7.2.2 Constructor & Destructor Documentation

#### 7.2.2.1 `__init__()`

```
def neutral_atom_imaging_simulation.Camera.CMOSCamera.__init__ (
    self,
    typing.Tuple[int,int] resolution,
    float dark_current_sampling_alpha = None,
    float dark_current_sampling_beta = None,
    float quantum_efficiency = None,
    float numerical_aperture = None,
    float physical_pixel_size = None,
    float magnification = None,
    float bias_clamp = None,
    float bias_stdev = None,
    float row_noise_stdev = None,
    float column_noise_scale = None,
    float flicker_noise_scale = None,
    float preampgain = None,
    float readout_stdev = None,
    float exposure_time = None,
    int binning = 1 )
```

#### 7.2.2.2 Constructor

Initializes all camera specific parameters and relays them to the library

## Parameters

<i>resolution</i>	Number of pixels per dimension
<i>dark_current_sampling_alpha</i>	Alpha value for sampling Gamma distribution for dark current
<i>dark_current_sampling_beta</i>	Beta value for sampling Gamma distribution for dark current
<i>quantum_efficiency</i>	Quantum efficiency of the camera for the corresponding wavelength [0.0,1.0]
<i>numerical_aperture</i>	Numerical apperture of the camera setup
<i>physical_pixel_size</i>	Physical pixel size of a single pixel (&#956m)
<i>magnification</i>	Magnification of the optical setup
<i>bias_clamp</i>	Bias clamp of the camera
<i>bias_stdev</i>	Standard deviation of the bias clamp
<i>row_noise_stdev</i>	Standard deviation of the row noise
<i>column_noise_scale</i>	Scale of the zero-mean Gumbel distribution for column noise
<i>flicker_noise_scale</i>	Scale of the zero-mean Gumbel distribution for flicker noise
<i>preampgain</i>	Preampgain of the camera
<i>readout_stdev</i>	Standard deviation of the final readout
<i>exposure_time</i>	Exposure time (s)
<i>binning</i>	Binning factor for the final image

## 7.2.3 Member Function Documentation

### 7.2.3.1 apply\_settings()

```
def neutral_atom_imaging_simulation.Camera.CMOSCamera.apply_settings (
    self )
```

Function for relaying any settings changes to the library.

## Returns

None

### 7.2.3.2 get\_image\_creation\_method()

```
def neutral_atom_imaging_simulation.Camera.CMOSCamera.get_image_creation_method (
    self )
```

Function for acquiring the function handle of the library that is used to generate images using this camera.

## Returns

The library function for generating images using this camera

Reimplemented from [neutral\\_atom\\_imaging\\_simulation.Camera.Camera](#).


The documentation for this class was generated from the following file:

- neutral\_atom\_imaging\_simulation/Camera.py

## 7.3 neutral\_atom\_imaging\_simulation.Camera.EMCCDCamera Class Reference

Use this camera if the generated images should look like they are taken by an electron multiplying charge-coupled device (EMCCD) camera.

Inheritance diagram for `neutral_atom_imaging_simulation.Camera.EMCCDCamera`:



## Public Member Functions

- `def __init__ (self, typing.Tuple[int, int] resolution, float dark_current_rate=None, float cic_chance=None, float quantum_efficiency=None, float numerical_aperture=None, float physical_pixel_size=None, float magnification=None, float bias_clamp=None, float preampgain=None, float scic_chance=None, float readout_stddev=None, int number_gain_reg=None, float p0=None, float exposure_time=None, int binning=1)`
- `def get_image_creation_method (self)`  
*Function for acquiring the function handle of the library that is used to generate images using this camera.*
- `def apply_settings (self)`  
*Function for relaying any settings changes to the library.*

## Public Attributes

- `dark_current_rate`
- `cic_chance`
- `quantum_efficiency`
- `numerical_aperture`
- `physical_pixel_size`
- `magnification`
- `bias_clamp`
- `preampgain`
- `scic_chance`
- `readout_stddev`
- `number_gain_reg`
- `p0`
- `exposure_time`
- `binning`
- `resolution`
- `zernike_coefficients`

### 7.3.1 Detailed Description

Use this camera if the generated images should look like they are taken by an electron multiplying charge-coupled device (EMCCD) camera.

### 7.3.2 Constructor & Destructor Documentation

### 7.3.2.1 `__init__()`

```
def neutral_atom_imaging_simulation.Camera.EMCCDCamera.__init__ (
    self,
    typing.Tuple[int,int] resolution,
    float dark_current_rate = None,
    float cic_chance = None,
    float quantum_efficiency = None,
    float numerical_aperture = None,
    float physical_pixel_size = None,
    float magnification = None,
    float bias_clamp = None,
    float preampgain = None,
    float scic_chance = None,
    float readout_stdev = None,
    int number_gain_reg = None,
    float p0 = None,
    float exposure_time = None,
    int binning = 1 )
```

### 7.3.2.2 Constructor

Initializes all camera specific parameters and relays them to the library

#### Parameters

<i>resolution</i>	Number of pixels per dimension
<i>dark_current_rate</i>	Rate of dark current (photons/s)
<i>cic_chance</i>	Chance for any given pixel to generate a clock-induced charge (photons)
<i>quantum_efficiency</i>	Quantum efficiency of the camera for the corresponding wavelength [0.0,1.0]
<i>numerical_aperture</i>	Numerical apperture of the camera setup
<i>physical_pixel_size</i>	Physical pixel size of a single pixel (&#956m)
<i>magnification</i>	Magnification of the optical setup
<i>bias_clamp</i>	Bias clamp of the camera
<i>preampgain</i>	Preampgain of the camera
<i>scic_chance</i>	Chance for a serial clock-induced charge to occure in any given gain register [0.0,1.0]
<i>readout_stdev</i>	Standard deviation of the final readout
<i>number_gain_reg</i>	Number of gain registers
<i>p0</i>	Chance for generating a secondary electron in any given gain register $(1 + p0)^{\text{number\_gain\_reg}} = \text{gain}$ [0.0,1.0]
<i>exposure_time</i>	Exposure time (s)
<i>binning</i>	Binning factor for the final image

## 7.3.3 Member Function Documentation

### 7.3.3.1 `apply_settings()`

```
def neutral_atom_imaging_simulation.Camera.EMCCDCamera.apply_settings (
    self )
```

Function for relaying any settings changes to the library.

**Returns**

None

**7.3.3.2 get\_image\_creation\_method()**

```
def neutral_atom_imaging_simulation.Camera.EMCCDCamera.get_image_creation_method (
    self )
```

Function for acquiring the function handle of the library that is used to generate images using this camera.

**Returns**

The library function for generating images using this camera

Reimplemented from [neutral\\_atom\\_imaging\\_simulation.Camera.Camera](#).


The documentation for this class was generated from the following file:

- [neutral\\_atom\\_imaging\\_simulation/Camera.py](#)

## 7.4 neutral\_atom\_imaging\_simulation.Experiment.Experiment Class Reference

Abstract experiment class.

Inheritance diagram for `neutral_atom_imaging_simulation.Experiment.Experiment`:

**Public Member Functions**

- def [get\\_atom\\_sites](#) (self)
- def **uses\_camera\_coords** (self)
- def [set\\_library](#) (self, ctypes.CDLL [library](#))  
*Function for setting the image generation library.*
- def [library](#) (self)  
*Function for getting the image generation library.*

**7.4.1 Detailed Description**

Abstract experiment class.

**7.4.2 Member Function Documentation****7.4.2.1 get\_atom\_sites()**

```
def neutral_atom_imaging_simulation.Experiment.Experiment.get_atom_sites (
    self )
```

Reimplemented in [neutral\\_atom\\_imaging\\_simulation.Experiment.TweezerArray](#).

### 7.4.2.2 library()

```
def neutral_atom_imaging_simulation.Experiment.Experiment.library (
    self )
```

Function for getting the image generation library.

#### Returns

The image generation C library

### 7.4.2.3 set\_library()

```
def neutral_atom_imaging_simulation.Experiment.Experiment.set_library (
    self,
    ctypes.CDLL library )
```

Function for setting the image generation library.

#### Parameters

<i>library</i>	The image generation C library
----------------	--------------------------------

#### Returns

None

The documentation for this class was generated from the following file:

- neutral\_atom\_imaging\_simulation/Experiment.py

## 7.5 neutral\_atom\_imaging\_simulation.ImageGenerator.ImageGenerator Class Reference

Main class for generating images.

### Public Member Functions

- def **\_\_init\_\_** (self)  
*Constructor Loads C library for later use.*
- def [get\\_library](#) (self)  
*Returns the loaded C library.*
- def [set\\_camera](#) (self, [Camera](#) camera)  
*Function for setting the camera.*
- def [set\\_experiment](#) (self, [Experiment](#) experiment)  
*Function for setting the experiment to be imaged.*
- def [create\\_image](#) (self)  
*Function to be called for generating an image.*
- def **read\_config\_file** (self, str path)

### 7.5.1 Detailed Description

Main class for generating images.

### 7.5.2 Member Function Documentation

### 7.5.2.1 create\_image()

```
def neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator.create_image (
    self )
```

Function to be called for generating an image.

#### Returns

Numpy array of generated image  
Numpy array of ground truths per atom site

### 7.5.2.2 get\_library()

```
def neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator.get_library (
    self )
```

Returns the loaded C library.

#### Returns

The image generation C library

### 7.5.2.3 set\_camera()

```
def neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator.set_camera (
    self,
    Camera camera )
```

Function for setting the camera.

#### Parameters

<i>camera</i>	The camera to be used for imaging
---------------	-----------------------------------

#### Returns

None

### 7.5.2.4 set\_experiment()

```
def neutral_atom_imaging_simulation.ImageGenerator.ImageGenerator.set_experiment (
    self,
    Experiment experiment )
```

Function for setting the experiment to be imaged.

#### Parameters

<i>experiment</i>	The experiment to be imaged
-------------------	-----------------------------

#### Returns

None

The documentation for this class was generated from the following file:


- neutral\_atom\_imaging\_simulation/ImageGenerator.py



## 7.6 neutral\_atom\_imaging\_simulation.Experiment.TweezerArray Class Reference

Use this if a tweezer array is to be simulated.

Inheritance diagram for neutral\_atom\_imaging\_simulation.Experiment.TweezerArray:



### Public Member Functions

- def `__init__` (self, stray\_light\_rate=None, imaging\_wavelength=None, scattering\_rate=None, survival\_probability=None, fill\_rate=None)  
*Constructor.*
- def `apply_settings` (self)  
*Function for relaying any settings changes to the library.*
- def `set_atom_sites_camera_space` (self, atom\_sites)  
*Function for setting the list of atom sites in normalized camera coordinates.*
- def `configure_atom_sites_camera_space` (self, tuple spacing, tuple count, tuple offset, float angle)  
*Function for configuring the layout of atom sites in normalized camera coordinates.*
- def `set_atom_sites_physical_space` (self, atom\_sites)  
*Function for setting the list of atom sites in physical object space.*
- def `configure_atom_sites_physical_space` (self, tuple spacing, tuple count, tuple offset, float angle)  
*Function for configuring the layout of atom sites in physical object space.*
- def `get_atom_sites` (self)  
*Function for acquiring the list of atom sites.*
- def `uses_camera_coords` (self)

### Public Attributes

- `stray_light_rate`
- `imaging_wavelength`
- `scattering_rate`
- `survival_probability`
- `fill_rate`
- `atom_sites`

#### 7.6.1 Detailed Description

Use this if a tweezer array is to be simulated.

#### 7.6.2 Constructor & Destructor Documentation

### 7.6.2.1 `__init__()`

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.__init__ (
    self,
    stray_light_rate = None,
    imaging_wavelength = None,
    scattering_rate = None,
    survival_probability = None,
    fill_rate = None )
```

Constructor.

#### Parameters

<i>stray_light_rate</i>	Rate of stray light (photons/s)
<i>imaging_wavelength</i>	The imaging wavelength (um)
<i>scattering_rate</i>	The average of photons emitted by an atom per second (photons/s)
<i>survival_probability</i>	The chance for an atom to survive being imaged [0.0,1.0]
<i>fill_rate</i>	The chance for an atom site to be filled

## 7.6.3 Member Function Documentation

### 7.6.3.1 `apply_settings()`

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.apply_settings (
    self )
```

Function for relaying any settings changes to the library.

#### Returns

None

### 7.6.3.2 `configure_atom_sites_camera_space()`

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.configure_atom_sites_camera_space
(
    self,
    tuple spacing,
    tuple count,
    tuple offset,
    float angle )
```

Function for configuring the layout of atom sites in normalized camera coordinates.

On average, there will be `column_count * row_count * fill_rate` atoms.

#### Parameters

<i>spacing</i>	The distance between two columns of atoms in normalized camera coordinates (x,y)
<i>count</i>	The number of columns (x,y)
<i>offset</i>	Distance from the image edge to the first atom site in normalized camera coordinates (x,y)
<i>angle</i>	The whole grid is rotated around the first atom site by this angle. Specify in radians

**Returns**

None

**7.6.3.3 configure\_atom\_sites\_physical\_space()**

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.configure_atom_sites_physical_↵
space (
    self,
    tuple spacing,
    tuple count,
    tuple offset,
    float angle )
```

Function for configuring the layout of atom sites in physical object space.  
On average, there will be `column_count * row_count * fill_rate` atoms.

**Parameters**

<i>spacing</i>	The distance between two columns of atoms in &#956m (x,y)
<i>count</i>	The number of columns (x,y)
<i>offset</i>	Distance from the image edge to the first atom site in &#956m (x,y)
<i>angle</i>	The whole grid is rotated around the first atom site by this angle. Specify in radians

**Returns**

None

**7.6.3.4 get\_atom\_sites()**

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.get_atom_sites (
    self )
```

Function for acquiring the list of atom sites.

**Returns**

The list of sites

Reimplemented from [neutral\\_atom\\_imaging\\_simulation.Experiment.Experiment](#).

**7.6.3.5 set\_atom\_sites\_camera\_space()**

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.set_atom_sites_camera_space (
    self,
    atom_sites )
```

Function for setting the list of atom sites in normalized camera coordinates.

**Parameters**

<i>atom_sites</i>	List of coordinate tuples where atom may reside, each coordinate value should be within [0.0,1.0]
-------------------	---

**Returns**

None

### 7.6.3.6 set\_atom\_sites\_physical\_space()

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.set_atom_sites_physical_space (
    self,
    atom_sites )
```

Function for setting the list of atom sites in physical object space.

#### Parameters

<i>atom_sites</i>	List of coordinate tuples where atom may reside, each coordinate value should be within [0.0,1.0]
-------------------	---

#### Returns

None

### 7.6.3.7 uses\_camera\_coords()

```
def neutral_atom_imaging_simulation.Experiment.TweezerArray.uses_camera_coords (
    self )
```

Reimplemented from [neutral\\_atom\\_imaging\\_simulation.Experiment.Experiment](#).

The documentation for this class was generated from the following file:

- neutral\_atom\_imaging\_simulation/Experiment.py